

Aeon7200

Anaesthesia Machine



- 7" color LCD screen, different parameters distinguished by colors for easy recognition
- Novel designed ventilator pneumatic system decreases the consumption of Oxygen
- Built-in "No Water Vapor Coagulation" sampling device ensures ZERO maintenance.
- Fast-release function by APL valve is comparable to high end models
- Advanced one button switch of manual /electronic provides convenient pneumatic /electronic controlling interchanges
- Optional AGSS device for more protection

Technical Specification

Physical & Environmental Specification

- Dimensions(H×L×W): 1386×839×628mm
- Weight: 85kg
- Storage temperature: -20~+55℃
- Working temperature: +10~+40℃
- Storage humidity: ≤93%

Electrical Specification

- Mains: 100-240V AC, 50/60Hz
- Battery: 24V DC, minimum 120 min
- Power consumption: ≤0.5A

Pneumatic Specification

- Gas supply: O₂, N₂O: 0.28-0.6Mpa
- Flowmeters: O₂: 0-1L/min; 1-10L/min
N₂O: 0-1L/min; 1-10L/min
- Gas system: O₂ deficiency alarm
Hypoxic guard system
O₂ flush: 25-75L/min
- Working mode: Closed, Semi-closed, Semi-open
- Driven mode: Pneumatically driven and electronically controlled
- Safety valve: ≤8kpa
- Operating mode: Man. / Vent.

Ventilator Specification

- Patient type: Adult, pediatric
- Ventilation modes: IPPV, Manual
- Setting: Touch key & navigator wheel knob
- Tidal volume: 0-1600ml
- Ventilation frequency: 4-100bpm
- I:E ratio: 4:1-1:8
- PEEP: 0-20 cmH₂O
- Inspiratory plateau: off, 5%-60%Ti
- High pressure: 21-80cmH₂O
- Low pressure: 0-20cmH₂O
- Monitoring display: Display: 7" color LCD screen
Pressure: P_{peak}, P_{plat}, P_{mean}, PEEP
Volumes: MV, V_t
Breath rate
Gas: FiO₂
Graph display: P-T, F-T waveforms
- Alarm: Pressure: P_{aw} high/low limit
Volumes: MV high/low limit
Gas: FiO₂ high/low limit
Airway pressure high
Power failure
Low O₂ supply pressure

Vaporizer

- Support 2 vaporizer (Selectatec compatible)
- Agent type: Halothane, Enflurane, Isoflurane, Sevoflurane available

Remark: Above configurations include standard and option. Please check price with your Aeonmed sales representative.

aeonmed
Beijing Aeonmed Co., Ltd.

Office (Headquarters):
Add: 11B2, Fengtai Science Park, (100070) Beijing, China
Tel: +86-10-8368 1616 Fax: +86-10-6371 8989
http://www.aeonmed.com E-mail: int@aeonmed.com

Information contained in this leaflet is correct at the date of publication. The policy of Beijing Aeonmed Co., Ltd. is one of the continued improvement to its products. Because of this policy Beijing Aeonmed Co., Ltd. reserves the right to make any change, which may affect the information in this leaflet without giving prior notices.

Print in July, 2012